Instruction Bulletin

Replaces MED10112110, 04/2011

ENGLISH

Iso-Gard™ IG2000 Series Remote Indicators for Line Isolation Monitors

Retain for future use.

Introduction

This bulletin provides installation, connection, and operation instructions for the Iso-Gard™ IG2000P and IG2000CBM remote indicators from Schneider Electric.

IG2000P

Product Description

The Iso-Gard™ IG2000P remote indicator from Schneider Electric provides remote indication of the visible and audible alarms from a line isolation monitor (LIM).

- Green LED—stays illuminated while the system is in normal condition
- Red hazard LED—illuminates when the Total Hazard Current exceeds the preset alarm level
- Audible hazard alarm—sounds when the Total Hazard Current exceeds the preset alarm level
- Mute button with yellow LED—silences the audible alarm on the remote indicator (local muting), or silences all audible alarms in the system (system muting)
- Test button—remotely performs a functional test of the LIM

The IG2000P remote indicator is available mounted in a one- or two-gang stainless steel faceplate for flush mounting into a panel or wall box with a two-inch minimum depth. The basic electrical connection is made using three wires. The wiring diagram on page 9 shows the connection to the Iso-Gard Series 6 (IG6) LIM.

Features




- Interfaces with the Iso-Gard LIM from Schneider Electric, as well as with LIMs from other manufacturers
- Does not interfere with medical equipment
- Uses low voltage wiring (12 Vdc or 12 Vac)
- Mounts to a standard, electrical gang box
- Connects by means of a screw terminal strip
- Uses long-life, LED displays
- Features an easy-to-clean, rugged, stainless steel and Lexan design

Configurations


The IG2000P Remote Indicator can be used in the following configurations:

- Wall-mounted
- Annunciator Stations
- Power Modules

Wall-Mounted

Description	Catalog Number	Mounting
Remote indicator mounted on a single-gang plate 	IG2000PG1	Customer-supplied, single-gang box
Remote indicator mounted on a two-gang plate 	IG2000PG2	Customer-supplied, two-gang box
Two remote indicators mounted on a two-gang plate 	IG2000P2G2	Customer-supplied, two-gang box

Annunciator Station

Description	Catalog Number▲	Mounting
Annunciator station 	SRAS()F	Back Box 2–8 remotes: SB120804 9–12 remotes: SB121204 13–16 remotes: SB121604 17+ remotes: quoted from factory

▲ To complete the catalog number, replace () with two digits (02–99) representing the number of IG2000P remote indicators to be installed onto one trim assembly.

Power Module

Description	Catalog Number ■	Mounting
Power receptacle module (208 or 240 V) 	SXRM()A1F SXRM()A2F	Back Box SB120804

■ To complete the catalog number, see the current version of catalog #4800CT1201, *Medical Isolated Power Panels*. Catalog numbers ending in 2F are equipped with an “In Use” light.

Safety Precautions

⚠ DANGER


HAZARD OF ELECTRIC SHOCK, EXPLOSION, OR ARC FLASH

- Apply appropriate personal protective equipment (PPE) and follow safe electrical work practices. See NFPA 70E or CSA Z462.
- This equipment must be installed and serviced only by qualified electrical personnel.
- Turn off all power supplying this equipment before working on or inside equipment.
- Always use a properly rated voltage sensing device to confirm power is off.


Failure to follow these instructions will result in death or serious injury.

Installation


IG2000PG1


IG2000PG2


IG2000P2G2


SRAS08F

SRAS12F

SRAS16F

Operation

The IG2000P remote indicator is controlled by a 12 Vdc signal from the LIM. As long as the system is in normal condition, the green LED stays illuminated.

When the LIM goes into an alarm condition, the hazard line is energized—the green LED turns off, the red LED illuminates, and the audible alarm sounds.

Pressing the Mute button mutes the audible alarm signal; the yellow LED illuminates to show that the audible alarm has been silenced. When the leakage current returns to an acceptable limit level, the red and yellow LEDs turn off and the green LED illuminates.


Pressing the Test button puts the LIM through a functional test cycle. The button must be pressed for approximately ten seconds to complete the test.

Technical Information

Item		Value
Operating voltage		12 Vdc or 12 Vac
Maximum current		40 mA
Operating class		Continuous operation
Ambient temperature	When operating	+32 °F to +122 °F (0 °C to +50 °C)
	When stored	-13 °F to +158 °F (-25 °C to +70 °C)
Connection		Screw terminal block
Maximum conductor size		16 AWG, stranded
Mounting		Screws
Weight	One-gang	0.25 lb. (0.11 kg)
	Two-gang	0.32 lb. (0.15 kg)

Connection to Iso-Gard Series 6 LIM

Wiring Diagram

Legend	Description
—	Connections required in systems using IG2000P
...	Additional connections required for "System Muting" option.
---	Additional connections required for "Test" option.
L1, L2	Connection to the isolated power system to be monitored.
12VDC CM	Common connection for Iso-Gard alarm indicator and test combinations, up to four (200 mA)
A, B	RS-485 interface, communication bus
RI1	+12 V connection for the test button of the Iso-Gard alarm indicator and test combination or voltage supply of IG2000P
K1/NC K1/COM K1/NO	N/C contact of alarm relay K1 Common contact of K1 N/O contact of K1
SAFE	System Normal signal for Iso-Gard alarm indicator and test combination that does not have bus compatibility
HAZARD	System Hazard signal for Iso-Gard alarm indicator and test combination that does not have bus compatibility
RI2	MUTE function for Iso-Gard alarm indicators that do not have bus compatibility
GND2, LIMGND	Separate connections to ground
TEST	Connection for testing the Iso-Gard LIM from any Iso-Gard alarm indicator that does not have bus compatibility
1S1, 1S2	Connection for current transformer CT1 ¹
Z1/M+, Z2/M-	Connection Z1/Z2 for over-temperature sensor or connection M+/M- for external measuring instrument for indication of Total Hazard Current
K2/COM K2/NC K2/NO	Common contact of alarm relay K2 N/C contact of K2 N/O contact of K2

¹ Conductors being monitored will be placed centrally through the current transformer.

Upgrading Legacy Remote Indicators

If you have this LIM	and this legacy remote alarm indicator,	you can upgrade to this remote alarm indicator.
 Iso-Gard Series 5		 IG2000PPM (includes IG2000P and panel mount bracket)
	 RA1 WM (wall-mounted)	 IG2000PG1 (designed for mounting in a single-gang electrical box)
	 IA1C	 IG2000PG2 (designed for mounting in a double-gang electrical box)
	 M5IAI	 IG2000PST0614 (designed for mounting in a legacy back box—530087BB)
	 M5IAI50	 IG2000PG4 (designed for mounting in a four-gang electrical box)

Connection to Iso-Gard Series D LIM



⚠ DANGER

HAZARD OF ELECTRIC SHOCK, EXPLOSION, OR ARC FLASH


- Apply appropriate personal protective equipment (PPE) and follow safe electrical work practices. See NFPA 70E.
- This equipment must be installed and serviced only by qualified electrical personnel.
- Turn off all power supplying this equipment before working on or inside equipment.
- Always use a properly rated voltage sensing device to confirm power is off.

Failure to follow these instructions will result in death or serious injury.

1. Make the wiring connections between the Series D connector plate and the IG2000P remote indicator as shown below. Use 16 ga., stranded wire.

Front View


Rear View

Connector Plate Terminal	Connects to Remote Terminal
1	2
2	1
3	3

IG2000CBM

Product Description

The Iso-Gard™ IG2000CBM remote indicator from Schneider Electric provides remote indication of the visible and audible alarms and digital mA reading from an Iso-Gard Series 6 (IG6) line isolation monitor (LIM).

- Green LED—stays illuminated while the system is in normal condition
- Red hazard LED—illuminates when the Total Hazard Current (THC) exceeds the preset alarm level
- Audible hazard alarm—sounds when the THC exceeds the preset alarm level
- Mute button with yellow LED—silences the audible alarm on the remote indicator (local muting), or silences all audible alarms in the system (system muting)
- Test button—remotely performs a functional test of the LIM
- Digital display—shows the THC in real time

The IG2000CBM remote indicator is available mounted in a two-gang stainless steel faceplate for flush mounting into a wall box with a two-inch minimum depth. The basic electrical connection is made using four wires. (See the wiring diagram on page 16.)

Features

- Interfaces with the Iso-Gard Series 6 LIM from Schneider Electric
- Does not interfere with medical equipment
- Uses low voltage wiring and RS-485 wiring
- Mounts to a standard, electrical two-gang box
- Connects by means of a screw terminal strip
- Uses long-life, LED displays
- Features an easy-to-clean, rugged, stainless steel and Lexan design

Configuration

The IG2000CBM remote indicator is available in a wall-mounted configuration.

Wall-Mounted

Description	Catalog Number	Mounting
Remote indicator mounted on a two-gang plate	IG2000CBM	Customer-supplied, two-gang box

Safety Precautions

⚠ DANGER


HAZARD OF ELECTRIC SHOCK, EXPLOSION, OR ARC FLASH

- Apply appropriate personal protective equipment (PPE) and follow safe electrical work practices. See NFPA 70E or CSA Z462.
- This equipment must be installed and serviced only by qualified electrical personnel.
- Turn off all power supplying this equipment before working on or inside equipment.
- Always use a properly rated voltage sensing device to confirm power is off.

Failure to follow these instructions will result in death or serious injury.

Installation

IG2000CBM

Operation

The IG2000CBM remote indicator is controlled by a 12 Vdc signal from the IG6 LIM. As long as the system is in normal condition, the green LED stays illuminated. Should the IG2000CBM lose communication with the IG6 LIM, the remote indicator will display 'ER' and the audible alarm will sound.

When the IG6 LIM goes into an alarm condition—the green LED turns off, the red LED illuminates, and the audible alarm sounds.

Pressing the Mute button mutes the audible alarm signal; the yellow LED illuminates to show that the audible alarm has been silenced. When the leakage current returns to an acceptable limit level, the red and yellow LEDs turn off and the green LED illuminates.

Pressing the Test button puts the IG6 LIM through a functional test cycle.


NOTE: The button must be pressed for approximately ten seconds to complete the test.

Technical Information

Item		Value
Operating voltage		12 Vdc
Maximum current		100 mA
Operating class		Continuous operation
Ambient temperature	When operating	+32 °F to +122 °F (0 °C to +50 °C)
	When stored	-13 °F to +158 °F (-25 °C to +70 °C)
Connection		Screw terminal block
Maximum conductor	Size	RS-485 twisted, shielded pair: 0.008–0.13 sq. in. (0.2–3.3 mm ²)
	Length	98.4 ft. (30 m)
Tightening torque		4.4–7.1 lb. in. (0.5–0.8 N•m)
Mounting		Screws
Weight		0.33 lb. (150 g)

Connection to Iso-Gard Series 6 LIM

Wiring Diagram

Legend	Description
—	Connections required in systems using IG2000CBM.
L1, L2	Connection to the isolated power system to be monitored
12VDC CM	Common connection for Iso-Gard alarm indicator and test combinations, up to two (200 mA)
A, B	RS-485 interface, communication bus
RI1	+12 V connection for the test button of the Iso-Gard alarm indicator and test combination or voltage supply of IG2000CBM
K1/NC	N/C contact of alarm relay K1
K1/COM	Common contact of K1
K1/NO	N/O contact of K1
SAFE	System Normal signal for Iso-Gard alarm indicator and test combination that does not have bus compatibility
HAZARD	System Hazard signal for Iso-Gard alarm indicator and test combination that does not have bus compatibility
RI2	MUTE function for Iso-Gard alarm indicators that do not have bus compatibility
GND2, LIMGND	Separate connections to ground
TEST	Connection for testing the Iso-Gard LIM from any Iso-Gard alarm indicator that does not have bus compatibility
1S1, 1S2	Connection for current transformer CT1 ¹
Z1/M+, Z2/M –	Connection Z1/Z2 for over-temperature sensor or connection M+/M- for external measuring instrument for indication of Total Hazard Current
K2/COM	Common contact of alarm relay K2
K2/NC	N/C contact of K2
K2/NO	N/O contact of K2

¹ Conductors being monitored will be placed centrally through the current transformer.

Connecting Multiple Remotes to a Single Iso-Gard Series 6 LIM


Up to two IG2000CBM remote indicators can be connected to a single IG6 LIM. IG2000CBM remote indicators use a special RS-485 protocol for communication. Each remote must be assigned a unique bus address for proper operation.

NOTE: The connected IG6 LIM will be assigned address 1. A single connected IG2000CBM will be assigned address 2. This is the factory default and will not require the procedure detailed in Step 2 below. An additionally connected remote indicator must be assigned address 3 and will require the procedure detailed in Step 2.

Step 1: Wiring Modifications

All remote indicators and the IG6 LIM are wired in series with each other via each device's "A" and "B" terminals. **Only** devices at the beginning and the end of the chain require the termination resistor. For the IG6 LIM, this internal resistor is activated via a switch on the device. For IG2000CBM remotes, a 120 Ω resistor must be wired across the remote's "A" and "B" terminals. The following figure shows a sample layout.

Sample Communication Bus Layout—Two Remote Indicators

Step 2: Changing the Address of the IG2000CBM

Follow this procedure to change the communication address of the IG2000CBM remote indicator.

SP

SP

Ad

02

Ad

03

Ar

01

1. Hold the TEST and MUTE buttons simultaneously for at least four seconds. "SP" will display on the upper and lower screens.
2. Push the TEST button until "Ad" is displayed in the top screen and two numbers (the factory default of "02") is displayed in the bottom screen. This number is the communication address of the IG2000CBM remote indicator.
3. Press the MUTE button. Press the TEST button until the desired address is displayed on the screen. Press the MUTE button to confirm this value.
4. Press the TEST button. "Ar" will display on the top screen, and two numbers (factory default "01") will display on the bottom screen. This is the address of the connected IG6 device. In most applications, this value will not need to be changed. If so, skip to step 5. Otherwise, press the MUTE button. Press the TEST button until the desired IG6 address appears, and press MUTE again to confirm this value.
5. Cycle power to the IG2000CBM to finish.

Upgrading Legacy Remote Indicators

If you have these LIM components	and this legacy remote alarm indicator,	you can upgrade to this remote alarm indicator.
 IG6M  IG6C  IG6BKTVM	 M5IAI  M5IAI50	 IG2000CBMST0614 (designed for mounting in a legacy back box-530087BB)  IG2000CBMG4 (designed for mounting in a four-gang electrical box)

Schneider Electric USA, Inc.
1415 S. Roselle Road
Palatine, IL 60067 USA
1-888-778-2733
www.schneider-electric.us

Electrical equipment should be installed, operated, serviced, and maintained only by qualified personnel. No responsibility is assumed by Schneider Electric for any consequences arising out of the use of this material.

© 2011–2014 Schneider Electric All Rights Reserved
All trademarks are owned by Schneider Electric Industries SAS or its affiliated companies.

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California Proposition 65 Warning—Nickel Compounds and Bisphenol A (BPA)

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⚠ WARNING: This product can expose you to chemicals including Nickel compounds, which are known to the State of California to cause cancer, and Bisphenol A (BPA), which is known to the State of California to cause birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov.

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All trademarks are the property of Schneider Electric SE, its subsidiaries, and affiliated companies.

Schneider Electric USA, Inc.
800 Federal Street
Andover, MA 01810 USA
888-778-2733
www.schneider-electric.us

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Toutes les marques commerciales sont la propriété de Schneider Electric Industries SAS ou de ses compagnies affiliées.
Schneider Electric México, S.A. de C.V.
Schneider Electric Canada, Inc.